



# Surgical Endoscopes All Medical Specialties


المناظير الجراحية الألمانية  
لكل التخصصات الطبية



الماسة العربية الزرقاء الطبية  
Arabian Blue Diamond



True vision difference.





**Arabian Blue Diamond**

Al Faisaliah Tower, Olaya District, Riyadh  
P.O Box 12234, Riyadh 3124 Kingdom of Saudi Arabia  
Tel: +966 (11) 490 3741  
Fax: +966 (11) 490 3752




### Advantages of flexible Emos endoskopos



Modern design



Made in Germany



Brilliant image representation

High-quality components „made in Germany“  
High depth of field and homogeneous illumination



Stable guiding element

Increasing lifetime due to innovative construction method



Robust shaft

Very good material resistance during the preparation process



Focusable

Optimal image setting through the ergonomic focusing ring



Standard eye piece

for camera application. Compatible with different camera TV coupler



Universal high guide / cable connection

Compatible with adapters from different manufacturers.

### Advantages of Emos optics



Sapphire negative & sapphire window

Solded sapphire negative and sapphire window  
advantages: no glue, no leakage and no scratches



Locking

Compatible with sheaths from different manufacturer  
advantages: easy combination with your existing equipment



Autoclavable

134°C / 273 °F / 2,3 bar / 21 min



Material

Complete made from stainless steel for high corrosion resistance



3-tube-system

Optic system is additional protected



Laser welding

Trough the laser welding absolutely sealed



Image quality

HD image display in real reproduction of the colors and tissue structures



Image quality

High-resolution, vignetting-free and distortion-free image quality in 4K

#### Arabian Blue Diamond

Al Faisaliah Tower, Olaya District, Riyadh  
P.O Box 12234, Riyadh 3124 Kingdom of Saudi Arabia  
Tel: +966 (11) 490 3741  
Fax: +966 (11) 490 3752



Click for  
more details

## ADVANTAGES OF FLEXIBLE EMOS ENDOSCOPES

*Good to know* ▶▶



### MODERN DESIGN

MADE IN GERMANY



### BRILLIANT IMAGE REPRESENTATION

High-quality components „made in Germany“  
High depth of field and homogeneous illumination



### STABLE GUIDING ELEMENT

Increasing lifetime due to innovative construction method



### ROBUST SHAFT

Very good material resistance during the preparation process



### FOKUSABLE

Optimal image setting through the ergonomic focusing ring



### STANDARD EYE PIECE for camera application

Compatible with different camera TV coupler



### UNIVERSAL LIGHT GUIDE / CABLE CONNECTION

Compatible with adapters from different manufacturers

#### Arabian Blue Diamond

Al Faisaliah Tower, Olaya District, Riyadh  
P.O Box 12234, Riyadh 3124 Kingdom of Saudi Arabia  
Tel: +966 (11) 490 3741  
Fax: +966 (11) 490 3752



Bronchoscopy

مناظير القصبات الهوائية



emos  
technology



Manufacturer

Origin

Product

Arabic Name

System Photo

# Flexibles- Bronchoskop 5.2/6.0

Flexible-Bronchoscope 5.2/6.0



Technische Daten  
Technical Data

Model	/ model	FBS-5225	FBS-6028
Einführschlauch	/ insertion tube	Ø 5,2mm	Ø 6,0mm
Arbeitskanal	/ working channel	Ø 2,3mm	Ø 2,8mm
Nutzlänge	/ working length	600mm	600mm
Blickwinkel	/ field of view	100°	100°
Abwinklung auf/ab	/ angulation up/down	180°/130°	180°/130°
Zubehör	/ accessories	Absaugventil, Dichtigkeitstester, Koffer	Absaugventil, Dichtigkeitstester, Koffer
		Suction valve, Leak tester, suitcase	Suction valve, Leak tester, suitcase
Artikel Nr.	/ item no.	FBS-5225C <span style="color: green;">■</span>	FBS-6028C <span style="color: green;">■</span>
		FBS-5225B <span style="color: blue;">■</span>	FBS-6028B <span style="color: blue;">■</span>

### Arabian Blue Diamond

Al Faisaliah Tower, Olaya District, Riyadh  
P.O Box 12234, Riyadh 3124 Kingdom of Saudi Arabia  
Tel: +966 (11) 490 3741  
Fax: +966 (11) 490 3752



Click for  
more details

Otolaryngology

مناظير الأنف والأذن والحنجرة



emos  
technology



▶ Sinuskope / Otoskope /  
Laryngoskope  
Sinusscopes / Otosscopes / Laryngoscopes



▶ Naso-Pharyngo-  
Laryngoskop 2.8/3.2/3.4  
Naso-Pharyngo-Laryngoscope 2.8/3.2/3.4



**Arabian Blue Diamond**

Al Faisaliah Tower, Olaya District, Riyadh  
P.O Box 12234, Riyadh 3124 Kingdom of Saudi Arabia  
Tel: +966 (11) 490 3741  
Fax: +966 (11) 490 3752



Urology

مناظير المسالك



emos  
technology



► **Semi-Flexibles-  
Ureterorenoskop**  
Semi-Flexible-Ureterorenoscope



► **Flexibles  
Uretero-Renoskop 8.4 Fr.**  
Flexible Uretero-Renoscope 8.4 Fr.



**Arabian Blue Diamond**

Al Faisaliah Tower, Olaya District, Riyadh  
P.O Box 12234, Riyadh 3124 Kingdom of Saudi Arabia  
Tel: +966 (11) 490 3741  
Fax: +966 (11) 490 3752



Nephroscope

مناظير الكلي



emos  
technology



# Semi-Flexibles Nephroskop

## Semi-Flexible Nephroscope



### Arabian Blue Diamond

Al Faisaliah Tower, Olaya District, Riyadh  
P.O Box 12234, Riyadh 3124 Kingdom of Saudi Arabia  
Tel: +966 (11) 490 3741  
Fax: +966 (11) 490 3752



Click for  
more details



Gynecology

مناظير طب النساء



emos  
technology



► Zystoskope / Hysteroskope  
Cystoscopes / Hysteroscopes



► Flexibles-  
Zystoskop 5.5  
Flexible-Cystoscope 5.5



**Arabian Blue Diamond**

Al Faisaliah Tower, Olaya District, Riyadh  
P.O Box 12234, Riyadh 3124 Kingdom of Saudi Arabia  
Tel: +966 (11) 490 3741  
Fax: +966 (11) 490 3752



Click for  
more details

Laparoscopy

مناظير البطن



emos  
technology



# Laparoskop 4K

Laparoscope 4K

NEW

4K  
ENDOSCOPY



### Arabian Blue Diamond

Al Faisaliah Tower, Olaya District, Riyadh  
P.O Box 12234, Riyadh 3124 Kingdom of Saudi Arabia  
Tel: +966 (11) 490 3741  
Fax: +966 (11) 490 3752



Click for  
more details

Arthroscopy

مناظير طب العظام والمفاصل



emos  
technology



▶ Arthroskop  
Arthroscope



**Arabian Blue Diamond**

Al Faisaliah Tower, Olaya District, Riyadh  
P.O Box 12234, Riyadh 3124 Kingdom of Saudi Arabia  
Tel: +966 (11) 490 3741  
Fax: +966 (11) 490 3752



Click for  
more details