

Video Intubation Laryngoscopy Emergency, Ambulance, ICU

التنظير الحنجري بالفيديو
للتنبيب (الطوارئ-الإسعاف-العناية المركزة)



الماسة العربية الزرقاء الطبية
Arabian Blue Diamond


A distinctive approach to life-saving..



TRACHWAY[®]
VIDEO INTUBATION SERIES

Video Intubation Laryngoscopy

التنظير الحنجري بالفيديو للتببيب



Arabian Blue Diamond

Al Faisaliah Tower, Olaya District, Riyadh
P.O Box 12234, Riyadh 3124 Kingdom of Saudi Arabia
Tel: +966 (11) 490 3741
Fax: +966 (11) 490 3752



Click for
more details

Trachway[®] Handle Set TVI-4000



60°

Adjustable angle

Thumb operation

60° upward

45° downward

45°



Video output port



Micro SD card slot



Mini USB Charging port



TVI-4056-B

Alumium carrying case

L 508 W 270 H 100 (mm)

Weight 1.5kg



▶ **Video Light Stylet**


Light Source	Surround LEDs Guiding red LED (for Ø 5.0)
Material	S.S. 304
Focus Range	5-50mm
Field of view	FOV (D): 83°

- TVI-4102 Ø 5.0 L 345mm
- TVI-4105 Ø 5.0 L 420mm
- TVI-4106 Ø 4.0 L 300mm
- TVI-4108 Ø 4.0 L 420mm



Malleability

Bend for difficult airway case and nasal intubation



Advantage

- Patented guiding red LED
- 304 stainless steel
- Submersible for disinfection
- Seamless seal and foggy free



Quick switch connector

Click-in design for quick switch among blade, stylet, and soft style systems



Four Sizes



TVI-4102 Ø 5.0mm L: 345mm
-General ETT (5.5 above)
-Reinforce ETT

TVI-4105 Ø 5.0mm L: 420mm
-35 Fr. DLT, Preform ETT

TVI-4106 Ø 4.0mm L: 300mm
-Pediatric ETT (4.5 above)

TVI-4108 Ø 4.0mm L: 420mm
-32 Fr. DLT, Preform ETT





► **Video Flexible Stylet**

Light Source	Surround LEDs
Material	Medical Grand PU
Focus Range	5-50mm
Field of view	FOV (D): 83°

- TVI-4301 Ø 5.0 L 650mm
- TVI-4302 Ø 4.0 L 750mm



Flexible Tube

Not easy to break



Plug in & out connector

Fast and efficient
when implement
tube exchange



Two Sizes

TVI-4301 OD 5.0mm Lenth 750mm
TVI-4302 OD 4.0mm Lenth 750mm

Arabian Blue Diamond

Al Faisaliah Tower, Olaya District, Riyadh
P.O Box 12234, Riyadh 3124 Kingdom of Saudi Arabia
Tel: +966 (11) 490 3741
Fax: +966 (11) 490 3752





Video Flexible Blade

Light Source	LEDs
Material	S.S. 316
Focus Range	5-50mm
Field of view	FOV (D):83°

• TVL-4220	MILL #0	Neonatal
• TVL-4221	MILL #1	Infant
• TVL-4232	MAC #2	Pediatric
• TVL-4233	MAC #3	Adult
• TVL-4234	MAC #4	Adult



Durability

- Submersible for disinfection
- Seamless seal
- Foggy free



Quick switch connector

Click-in design for quick switch among blade, style, and soft style systems



Advantage

- Conventional blade design
- Same skill with non direct view
- 316 stainless steel



Five Sizes

TVL-4220 Mill#0	Neonatal
-W: 11.5mm L: 75mm	
TVL-4221 Mill#1	Infant
-W: 11.5mm L: 102mm	
TVL-4232 MAC#2	Pediatric
-W: 12mm L: 100mm	
TVL-4233 MAC#3	Adult
-W: 15mm L: 130mm	
TVL-4234 MAC#4	Adult
-W: 15mm L: 135mm	

Arabian Blue Diamond

Al Faisaliah Tower, Olaya District, Riyadh
P.O Box 12234, Riyadh 3124 Kingdom of Saudi Arabia
Tel: +966 (11) 490 3741
Fax: +966 (11) 490 3752

