



Surgical Endoscopes All Medical Specialties

المناظير الجراحية الألمانية
لكل التخصصات الطبية




ABD

إيه بي دي ميد إيست
ABD Mideast



True vision difference.





ABD MIDEAST GENERAL TRADING

Level 41, Emirates Towers Sheikh Zayed Road. Dubai, UAE PO Box 31303
Tel: +971 4 319 7867
Fax: +971 4 319 7897



Advantages of flexible Emos endoskopes



Modern design



Stable guiding element

Increasing lifetime due to innovative construction method




Made in Germany



Robust shaft

Very good material resistance during the preparation process



Brilliant image representation

High-quality components „made in Germany“
High depth of field and homogeneous illumination



Focusable

Optimal image setting through the ergonomic focusing ring



Standard eye piece

for camera application. Compatible with different camera TV coupler



Universal high guide / cable connection

Compatible with adapters from different manufacturers.

Advantages of Emos optics



Sapphire negative & sapphire window

Solded sapphire negative and sapphire window
advantages: no glue, no leakage and no scratches



Locking

Compatible with sheaths from different manufacturer
advantages: easy combination with your existing equipment



Autoclavable

134°C / 273°F / 2,3 bar / 21 min



Material

Complete made from stainless steel for high corrosion resistance



3-tube-system

Optic system is additional protected



Laser welding

Trough the laser welding absolutely sealed



Image quality

HD image display in real reproduction of the colors and tissue structures



Image quality

High-resolution, vignetting-free and distortion-free image quality in 4K

ADVANTAGES OF FLEXIBLE EMOS ENDOSCOPES

Good to know



MODERN DESIGN

MADE IN GERMANY



BRILLIANT IMAGE REPRESENTATION

High-quality components „made in Germany“
High depth of field and homogeneous illumination



STABLE GUIDING ELEMENT

Increasing lifetime due to innovative construction method



ROBUST SHAFT

Very good material resistance during the preparation process



FOKUSABLE

Optimal image setting through the ergonomic focusing ring



STANDARD EYE PIECE for camera application

Compatible with different camera TV coupler



UNIVERSAL LIGHT GUIDE / CABLE CONNECTION

Compatible with adapters from different manufacturers

ABD MIDEAST GENERAL TRADING

Level 41, Emirates Towers Sheikh Zayed Road. Dubai, UAE PO Box 31303

Tel: +971 4 319 7867

Fax: +971 4 319 7897



Bronchoscopy

مناظير القصبات الهوائية



emos
technology



Manufacturer

Origin

Product

Arabic Name

System Photo

Flexibles- Bronchoskop 5.2/6.0

Flexible-Bronchoscope 5.2/6.0



Technische Daten
Technical Data

Model	/ model	FBS-5223	FBS-6028
Einführschlauch	/ insertion tube	Ø 5,2mm	Ø 6,0mm
Arbeitskanal	/ working channel	Ø 2,3mm	Ø 2,8mm
Nutzlänge	/ working length	600mm	600mm
Blickwinkel	/ field of view	100°	100°
Abwinklung auf/ab	/ angulation up/down	180°/130°	180°/130°
Zubehör	/ accessories	Absaugventil, Dichtigkeitstester, Koffer	Absaugventil, Dichtigkeitstester, Koffer
		Suction valve, Leak tester, suitcase	Suction valve, Leak tester, suitcase
Artikel Nr.	/ item no.	FBS-5223C ■	FBS-6028C ■
		FBS-5223B ■	FBS-6028B ■

ABD MIDEAST GENERAL TRADING

Level 41, Emirates Towers Sheikh Zayed Road. Dubai, UAE PO Box 31303
Tel: +971 4 319 7867
Fax: +971 4 319 7897



Otolaryngology

مناظير الأنف والأذن والحنجرة



emos
technology



▶ Sinuskope / Otoskope /
Laryngoskope
Sinusscopes / Otosscopes / Laryngoscopes



▶ Naso-Pharyngo-
Laryngoskop 2.8/3.2/3.4
Naso-Pharyngo-Laryngoscope 2.8/3.2/3.4



ABD MIDEAST GENERAL TRADING

Level 41, Emirates Towers Sheikh Zayed Road. Dubai, UAE PO Box 31303
Tel: +971 4 319 7867
Fax: +971 4 319 7897





ABD

إيه بي دي ميد إيست
ABD Mideast

2024

New
Product

Urology

مناظير المسالك



emos
technology



▶ Semi-Flexibles-
Ureterorenoskop
Semi-Flexible-Ureterorenoscope



▶ Flexibles
Uretero-Renoskop 8.4 Fr.
Flexible Uretero-Renoscope 8.4 Fr.



ABD MIDEAST GENERAL TRADING

Level 41, Emirates Towers Sheikh Zayed Road. Dubai, UAE PO Box 31303
Tel: +971 4 319 7867
Fax: +971 4 319 7897

E-mail: info@abd-mideast.com

www.abd-mideast.com



Click for Links

Nephroscope

مناظير الكلي



emos
technology



Semi-Flexibles Nephroskop

Semi-Flexible Nephroscope



ABD MIDEAST GENERAL TRADING

Level 41, Emirates Towers Sheikh Zayed Road. Dubai, UAE PO Box 31303
Tel: +971 4 319 7867
Fax: +971 4 319 7897

E-mail: info@abd-mideast.com

www.abd-mideast.com



Click for Links

Gynecology

مناظير طب النساء



emos
technology



► Zystoskope / Hysteroskope
Cystoscopes / Hysteroscopes



► Flexibles-
Zystoskop 5.5
Flexible-Cystoscope 5.5



ABD MIDEAST GENERAL TRADING

Level 41, Emirates Towers Sheikh Zayed Road. Dubai, UAE PO Box 31303
Tel: +971 4 319 7867
Fax: +971 4 319 7897

E-mail: info@abd-mideast.com

www.abd-mideast.com



Click for Links

Laparoscopy

مناظير البطن



emos
technology



Laparoskop 4K

Laparoscope 4K

NEW

4K
ENDOSCOPY



ABD MIDEAST GENERAL TRADING

Level 41, Emirates Towers Sheikh Zayed Road. Dubai, UAE PO Box 31303
Tel: +971 4 319 7867
Fax: +971 4 319 7897

E-mail: info@abd-mideast.com

www.abd-mideast.com



Click for Links

Arthroscopy

مناظير طب العظام والمفاصل



emos
technology



▶ Arthroskop
Arthroscope



ABD MIDEAST GENERAL TRADING

Level 41, Emirates Towers Sheikh Zayed Road. Dubai, UAE PO Box 31303
Tel: +971 4 319 7867
Fax: +971 4 319 7897

E-mail: info@abd-mideast.com

www.abd-mideast.com



Click for Links