

Video Intubation Laryngoscopy Emergency, Ambulance, ICU

التنظير الحنجري بالفيديو
للتنبيب (الطوارئ-الإسعاف-العناية المركزة)



ABD

إيه بي دي ميد إيست
ABD Mideast

A distinctive approach to life-saving..





ABD

إيه بي دي ميد إيست
ABD Mideast


2024

New
Product

TRACHWAY®
VIDEO INTUBATION SERIES

Video Intubation Laryngoscopy

التنظير الحنجري بالفيديو للتببيب



Video Laryngoscope

ABD MIDEAST GENERAL TRADING

Level 41, Emirates Towers Sheikh Zayed Road. Dubai, UAE PO Box 31303

Tel: +971 4 319 7867

Fax: +971 4 319 7897


E-mail: info@abd-mideast.com

www.abd-mideast.com



Click for Links

Trachway® Handle Set
TVI-4000



Adjustable angle
Thumb operation
60° upward
45° downward



- Video output port
- Micro SD card slot
- Mini USB Charging port



TVI-4056-B Aluminum carrying case
L 508 W 270 H 100 (mm)
Weight 1.5kg



ABD

إيه بي دي ميد إيست
ABD Mideast

2024

New Product

TRACHWAY®
VIDEO INTUBATION SERIES

Video Malleable Stylet



▶ Video Light Stylet


Light Source	Surround LEDs Guiding red LED (for Ø 5.0)
Material	S.S. 304
Focus Range	5-50mm
Field of view	FOV (D): 83°

- TVI-4102 Ø 5.0 L 345mm
- TVI-4105 Ø 5.0 L 420mm
- TVI-4106 Ø 4.0 L 300mm
- TVI-4108 Ø 4.0 L 420mm



Malleability

Bend for difficult airway case and nasal intubation



Advantage

- Patented guiding red LED
- 304 stainless steel
- Submersible for disinfection
- Seamless seal and foggy free



Quick switch connector

Click-in design for quick switch among blade, stylet, and soft style systems



Four Sizes

TVI-4102 Ø 5.0mm L: 345mm
-General ETT (5.5 above)
-Reinforce ETT

TVI-4105 Ø 5.0mm L: 420mm
-35 Fr. DLT, Preform ETT

TVI-4106 Ø 4.0mm L: 300mm
-Pediatric ETT (4.5 above)

TVI-4108 Ø 4.0mm L: 420mm
-32 Fr. DLT, Preform ETT

ABD MIDEAST GENERAL TRADING

Level 41, Emirates Towers Sheikh Zayed Road. Dubai, UAE PO Box 31303

Tel: +971 4 319 7867

Fax: +971 4 319 7897

E-mail: info@abd-mideast.com

www.abd-mideast.com



Click for Links



ABD


إيه بي دي ميديا إيست
ABD Mideast

2024

New Product

TRACHWAY[®]
VIDEO INTUBATION SERIES

Video Flexible Stylet



► **Video Flexible Stylet**

Light Source	Surround LEDs
Material	Medical Grand PU
Focus Range	5-50mm
Field of view	FOV (D): 83°

- TVI-4301 Ø 5.0 L 650mm
- TVI-4302 Ø 4.0 L 750mm



Flexible Tube

Not easy to break



Plug in & out connector

Fast and efficient when implement tube exchange



Two Sizes

TVI-4301 OD 5.0mm Lenth 750mm
TVI-4302 OD 4.0mm Lenth 750mm

ABD MIDEAST GENERAL TRADING

Level 41, Emirates Towers Sheikh Zayed Road. Dubai, UAE PO Box 31303

Tel: +971 4 319 7867

Fax: +971 4 319 7897

E-mail: info@abd-mideast.com

www.abd-mideast.com



Click for Links



ABD

إيه بي دي ميد إيست
ABD Mideast

2024

New Product

TRACHWAY®
VIDEO INTUBATION SERIES

Video Rigid Laryngeal Blades



Video Flexible Blade

Light Source	LEDs
Material	S.S. 316
Focus Range	5-50mm
Field of view	FOV (D):83°

- TVL-4220 MILL #0 Neonatal
- TVL-4221 MILL #1 Infant
- TVL-4232 MAC #2 Pediatric
- TVL-4233 MAC #3 Adult
- TVL-4234 MAC #4 Adult



Durability

- Submersible for disinfection
- Seamless seal
- Foggy free



Quick switch connector

Click-in design for quick switch among blade, style, and soft style systems



Advantage

- Conventional blade design
- Same skill with non direct view
- 316 stainless steel



Five Sizes

TVL-4220 Mill#0	Neonatal
-W: 11.5mm L: 75mm	
TVL-4221 Mill#1	Infant
-W: 11.5mm L: 102mm	
TVL-4232 MAC#2	Pediatric
-W: 12mm L: 100mm	
TVL-4233 MAC#3	Adult
-W: 15mm L: 130mm	
TVL-4234 MAC#4	Adult
-W: 15mm L: 135mm	

ABD MIDEAST GENERAL TRADING

Level 41, Emirates Towers Sheikh Zayed Road. Dubai, UAE PO Box 31303
Tel: +971 4 319 7867
Fax: +971 4 319 7897

E-mail: info@abd-mideast.com

www.abd-mideast.com



Click for Links