

KAZAMA ENT TOTAL CATALOGUE





ENT treatment total unit
KAU-3000 Grand Type

KAU-3000

[HARMONY]

FEATURES

This product is comprehensive unit for ENT diagnosis and treatment with new and high technology functions, especially focused on more convenient examination and the best sanitary environment installing UV Sterilization system (option) for endoscope and other instruments as well. Anti-fog and endoscope holders are designed ergonomically to reduce the doctors' fatigue and the heating device (option) can be equipped to maintain always same temperature of instrument tray.

In addition, the visual system is mounted under the unit table to prevent obstruction of treatment causing from the wires.


Unit table is made of sustainable and durable artificial marble and special synthetic resin not only for luxurious looking, but also to prevent any kind of contamination from drug and other chemicals.

SPECIFICATIONS

Compressor pump	250W, Two head motor air pressure 7kgf/cm ²
Vacuum pump(main)	Negative pressure 680mmHg
Vacuum pump(sub)	Negative pressure 650mmHg
Main suction bottle	3000cc
Sub suction bottle	1000cc
Clean system	Air filter & Regulator
Power Consumption	600~1400W
Voltage	220V
Frequency	50~60Hz
Dimension	1950(W) x 620(L) x 740(H)
Weight	185kg
Anti-Fog device	400W



DETAIL VIEW



ENT treatment total unit
SPECIFICATION & STANDARD SET

2

NET-1500A,B,C

This is a patient chair, one of essential items of treatment, not only used in otorhinolaryngology hospital, but also the other medical sections such as pediatrics, family medicine and so on. It is designed to be adjustable to the heights and directions convenient for both of doctor's most effective treatment and patient's comfort.

FEATURES

- Hydraulic auto brake system(rotation)
- It can be turned more than 360°(unlimited)
- Water-proof foot switch : Protect from electric shock and erode problems
- Developed slide cylinder the first in Korea and it elevate durability more than 10 times

SPECIFICATIONS

Reclining angle	95°
Base Diameter	580mm
Power consumption	545W
Voltage	220V
Weight	95kg
Up-Down stroke	520mm(min) ~ 820mm(max)
Rotation angle	360°
Frequency	50 / 60 Hz
Dimension	620(W) x 650(L) x 1340(H) S/T300
Head rest	Manual type
Auto head rest	Option

SPECIFICATIONS

Foot paddle	1
Operating manual	1



ENT Chair

NET-1500 A,B,C



KAZAMA ENT

223-612, Seoknam-dong, Seo-gu, Incheon, KOREA, 404-220
Tel : 82-32-581-1900 Fax : 82-32-582-1309
<http://www.kazamaent.co.kr> <http://kazamaent.koreasme.com>