



TOTAL CATALOGUE

Our mission is to improve the health
and quality of life of those in need.

www.kazamaent.co.kr



KAU-3000[MINI]

FEATURES

This is a MINI Type product with new function of the otolaryngology examination and treatment unit.

SPECIFICATIONS

Compressor pump	70W 2~2.5kgf/cm ²
Vacuum pump(main)	Negative pressure 650mmHg
Main suction bottle	3000cc x 1
Clean system	Air filter & Regulator
Power Consumption	300W
Voltage	220V
Dimension(mm)	500(W) x 600(D) x 940(H)
Weight	30kg
Anti-fog device	250W

STANDARD ACC'S

Straight spray	2
Suction	1
Gauze container	2
Waste receptacle	2
Medicine bottle	2
instrument tray	1

OPTIONS

Att Arm Set	1
CCD camera & Light source	1
Doctor's stool	1
ENT chair	1



KAU-3000[MINI]



NET-1500[D,E]

It's an essential patient chair to treat of ear, nose and throat and it has several suitable functions depends on ears which located in both side, discrepancy in height and various position for treatment.

NET-1500[D] SPECIFICATIONS

Power Source	AC 220V, 50/60Hz
Power Consumption	300VA
Rotation Angle	340°
Return Position	Auto
Up-Down Stroke	250(480~730mm)
Reclinling Angle	90°
Weight	Approx. 75kg
Dimension(mm)	600(W) x 730(D) x 1080(H)
Option	Headrest



NET-1500[D]

NET-1500[E] SPECIFICATIONS

Power Source	AC 220V, 50/60Hz
Power Consumption	300VA
Rotation Angle	340°
Return Position	Auto
Up-Down Stroke	250(480~730mm)
Reclinling Angle	45°
Weight	Appro, 75kg
Dimension(mm)	600(W) x 730(D) x 1080(H)
Option	Headrest



NET-1500[E]



KAZAMA ENT

451, Hwahap-ro Beobwon-eup, Paju-si, Gyeonggi-do, KOREA (Zip code : 10807)

TEL 031-958-1902 | FAX 031-958-1302

www.kazamaent.co.kr www.kazamaent.com