

# MICRO-100(J)

## FEATURES

- With new optic fiber, can make more clear image
- It can be mounted at ENT basic unit
- CCD Camera can be attachable in order to see with the monitor(optional)

## SPECIFICATIONS

Focus length of objective	1F=250mm
Magnification	X2, X4, X7
Working distance	190mm
Bulb	12V/100W
Stand type or Unit attach type	
Parameter of Binocular tube	
- Inclined binocular tubes	45 ▲
- Diopter adjustment range	± 6D
- Pupil distance range	50mm / 70mm
- Height of eyecups	18mm
Parameter of illumination	
- 6 + 0 Coaxial illumination	
- Maximum intensity	30000Lx or up
Parameter of position adjustment	
- Maximum stretching radius of arm	870mm
- Vertical movement range	
(from floor to front surface of big objective)	700 / 1100mm
- Range of fine focusing adjustment	30mm
Parameter of Electric	
- Input voltage	AC 110/220V 50~60Hz
- Input power	120VA



MICROSCOPE UNIT

32